Schaper
HEALTHCARE

DE EN FR NL DK SE

OHS Knieschiene PCL

OHS Knee Brace PCL


Genouillère OHS PCL

OHS Kniebrace PCL

OHS Knäskinne PCL



OHS Knäskena PCL

Model: SC32I

DK

SE

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www.schaper-healthcare.com

Version V01/2022/01

Size	XS	S	M	L	XL	XXL
Thigh circumference (cm)	33-39	39-45	45-51	51-57	57-63	63-67
Side	Left/Right					

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Aluminium, Polyamid (PA), PUR-Schaumstoff, Polyethylen (PE), Polyester (PES)

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Intended use:

Knee joint orthosis for stabilising and guiding the knee joint with the option of restricting the range of movement (extension and flexion) of the joint.

Indications:

For treatment of the posterior cruciate ligament (PCL) after ligament ruptures and injuries
Ruptures and injuries of the collateral ligaments
Meniscus injuries
Postoperative treatment

Contraindications:
Currently unknown.**Adjusting the extension / flexion limit:**

You can adjust the desired range of movement at the joints of the brace. To do this, open the joint cover by loosening the two screws using the hex key (1). Take out the limit pins and insert them at the appropriate degree setting (2). Close the joint cover completely and tighten the screws again. Flexion limit: 0°, 20°, 45°, 60°, 90°
Extension limit: 0°, 5°, 20°, 30°, 45°

Application instruction:

If you can, sit on a chair and bend your leg slightly to make it easier to put the brace on. Open the back and front straps of the brace on one side. Place the brace around the leg (3). To do this, pass your foot through the frame and pull the brace up to the knee. Position the joints of the brace on the side of the knee so that the centre of the joint is at the level of the upper edge of the kneecap. The pads of the straps can be freely positioned using the Velcro fastener. When closing the straps, position the pads so that they are in the middle of the leg. First close the front top lower-leg strap and then the front bottom lower-leg strap (4+5). Then close the back bottom thigh strap and then the back top thigh strap (6). Finally, close the back lower-leg strap and the front thigh strap (7+8).

Application note / Caution:

- Users are generally advised not to wear the knee brace directly on the skin.
- In the case of skin problems (skin diseases, venous diseases, hypersensitivity to pressure, etc.), the knee brace should not be worn directly on the skin. For example, an anti-thrombosis stocking can be worn under the knee brace. In this case, consult your attending physician.
- The frame of the knee brace should not be too heavily deformed, otherwise the function of the joints may be impaired.
- When using the knee brace, the boundary pins must always be placed in the joint and must never be completely removed.
- In order to prevent the knee joint from spreading, a limit pin must always be placed in the degree setting 0° extension, unless your doctor specifies a different setting.
- Do not use the product on open wounds. In the following cases you should consult your doctor before use: Skin diseases, lymphatic drainage disorders, sensory disorders, circulatory disorders, unclear soft tissue swelling.
- This product is intended for use on a single patient only.
- If necessary, the help of a second person can be beneficial for the correct application of the product.
- Only by choosing the correct size and variant, as well as the correct application, an optimal function can be achieved.
- The product should no longer be used if it is damaged or heavily worn, or if existing joints and fasteners no longer function properly.
- Make sure that no fabric parts are damaged by existing Velcro fasteners.
- Do not put on the product too tightly, as this may affect the blood circulation.
- If symptoms persist or worsen, please consult your medical doctor.
- Please inform the manufacturer and the competent authority in your country about serious incidents with this product.
- If you are unclear or have questions about the product and its handling, please consult your doctor or qualified personnel.
- Do not make any changes to the product.
- Simultaneous use with other products only after consultation with your doctor.

Care instructions:

The textile parts of the product can be cleaned by hand with water and a mild detergent. The plastic and metal parts of the product can be cleaned with a damp cloth. Close existing Velcro fasteners before washing to avoid damaging other laundry items. Allow to air dry. Avoid direct heat sources and do not machine dry or iron. Make sure the product is completely dry before using it again. Do not dry clean and do not use bleach, solvents, or fabric softeners. Fatty or acidic agents, ointments or lotions can damage the material.

Material composition:

Aluminium, Polyamide (PA), PUR foam, Polyethylene (PE), Polyester (PES)

Storage instructions:

Store the product in a dry place and protect it from heat, direct sunlight and moisture.

Disposal note:

Disposal via household waste is possible.

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